



## IntelliPlex™ 10X Wash Buffer

**REF** 80210 100 mL

**IVD** For In Vitro Diagnostic Use



### PlexBio Co., Ltd.

6F-1, No. 351, Yangguang St., Neihu Dist.  
Taipei City, 11491, Taiwan  
TEL: + 886-2-2627-5878  
<http://www.plexbio.com>

### European Authorized Representative:



Medical Device Safety Services GmbH (MDSS)  
Schiffgraben 41, 30175, Hannover, Germany



**IMPORTANT:**  
Read the instructions carefully prior to use

### 1. INTENDED USE

The IntelliPlex™ 10X Wash Buffer enables users to easily prepare consumptive wash buffer for IntelliPlex™ 1000  $\pi$ Code Processor.

### 2. WARNINGS AND PRECAUTIONS

- For in-vitro diagnostic use.
- This wash buffer should be used by qualified laboratory personnel only.
- Do not use the buffer past its expiration date.
- Store the buffer at 2-8°C according to the product label and instructions.
  - General laboratory precautions should be taken: Do not pipette by mouth.
  - Wear protective clothing (e.g., disposable powderless gloves and laboratory coats) and eye protection.
  - Do not eat, drink or smoke in the laboratory.
  - Wash hands thoroughly after handling samples and reagents.
- Material Safety Data Sheets (MSDS) are available upon request from PlexBio Customer Service.

### 3. CONTENTS

- Phosphate buffered saline
- Tween-20 1 %
- Sodium azide <0.1 %

### 4. MATERIALS AND EQUIPMENT REQUIRED BUT NOT SUPPLIED

- Disposable gloves, powder-less
- ddH<sub>2</sub>O for dilution of 10X Wash Buffer
- U Tray (PlexBio; Cat. No. 80023)
- IntelliPlex 1000  $\pi$ Code Processor (PlexBio; Cat. No. 80033)

### 5. STORAGE, STABILITY AND TRANSPORTATION

The IntelliPlex 10X Wash Buffer should be kept at 2~8 °C. The shipping temperature for 10X Wash Buffer is under the room temperature. Do not use the IntelliPlex 10X Wash Buffer after the expiration date. The buffer is guaranteed up to the expiration date on the label if handled and stored under the recommended conditions.

### 6. INSTRUMENT AND BUFFER PREPARATION













Refer to the user manual of IntelliPlex 1000  $\pi$ Code Processor for complete operation instructions.

#### Buffer Preparation

Wash buffer purchased from PlexBio is 10X concentrated. Dilution is required before use. Transfer 100mL of the 10X Wash Buffer to the IntelliPlex 1000  $\pi$ Code Processor 1L Wash Buffer bottle and add 900 ml ddH<sub>2</sub>O. Mix by swirling.

**NOTE:** The prepared 1X Wash Buffer can be used for up to one week.


## 7. SYMBOLS

Symbol	Explanation	Symbol	Explanation
	In-vitro diagnostic use		Catalog number
	Serial number		Consult instructions for use
	Manufacturer		Use by Date
	Temperature limitation		Caution
	Manufacture Date		Batch number
	European Union Conformity		European Authorized Representative

### **Notice to User**

The use of this product and the associated PlexBio instrumentation is covered by one or more issued (US10302640B2, US10436778B2, US10436776B2, US9063044B2, US10019815B2) and pending US and foreign patents owned by PlexBio Co., Ltd. The purchase of this product includes nontransferable rights to use this amount of the product to practice the methods described therein. No general patent or other license of any kind other than this specific right of use from purchase is granted. Further information on purchasing licenses for other applications can be obtained from PlexBio Co., Ltd. 6F-1, No. 351, Yangguang St., Neihu District, Taipei City 11491, Taiwan.

### **Trademarks**

 ,  $\pi$ Code™ , DeXipher™, IntelliPlex™, DigiPlex™,  $\pi$ Code™ , PlexBio™, are trademarks or registered trademarks of PlexBio Co., Ltd. All other product names and any registered and unregistered trademarks mentioned are used for identification purposes only and remain the exclusive property of their respective owners.

### **Copyright**

©2020 PlexBio Co., Ltd. All rights reserved.