

Long-term Shutdown

Software version 1.1.8.1429



Follow this protocol if you will not be using the PlexBio 100 for more than two weeks.


Required item

- Intelliplex Calibration Kit (Ref. 80035-R)


PlexBio 100

Long-term Shutdown_A-3_V01

- Turn off the **PlexBio 100** and unplug all cables



- Please note: Calibration should be performed before the instrument is used again after long-term shutdown
(see **A-2_Performing Calibration** for details)



Thank You

PlexBio

6F-1, No. 351, Yangguang St.
Neihu District, Taipei City 11491, Taiwan
+866-2-2627-5878 | PlexBio.com
service@plexbio.com