

SAFETY DATA SHEET

Ver. 03 08/24/2020

SECT	ION 1: CHEMICAL PRODUCT AND COMPANY IDENTIFICATION
1.1. Product Identifier:	
Product Name	IntelliPlex RET/NTRK1 Rearrangement Kit
Product Code	82025
1.2. Relevant identified use	s of the substance or mixture and uses advised against
Identified uses	For professional users only.
1.3. Details of the supplier	of the safety data sheet
Manufacturer/Supplier:	PlexBio Co., Ltd.
	6F-1, No. 351 Yangguang Street, Neihu District, Taipei 114 Taiwan, R.O.C.
	Tel: +886-2-2627-5878 Fax: +886-2-2627-5979
1.4 Emergency telephone n	umber
Emergency Phone #	Tel: +886-2-2627-5878
Email	service@plexbio.com

SECTION 2: HAZARD IDENTIFICATION

2.1 Classification of the substance or mixture

The product is a kit consisting of individual ingredients listed in section 3. Not a hazardous substance or mixture according to Regulation (EC) No. 1272/2008 or 2012 OSHA Hazard Communication Standard (29 CFR 1910.1200). Not a dangerous substance or mixture according to the Globally Harmonized System (GHS).

2.2 GHS Label elements, including precautionary statements:

Not a hazardous substance or mixture.

2.3 Hazards not otherwise classified (HNOC) or not covered by GHS

None identified. This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent and very bioaccumulative (vPvB) at levels of 0.1% or higher.

SECTION 3: COMPOSITION AND INFORMATION ON INGREDIENTS			
Kit Contents	Known Hazardous Components	CAS-No.	Concentration
RET/NTRK1 KIT RT-PCR Buffer	Not Hazardous	-	-
RET/NTRK1 KIT RT-PCR PrimerMix	Not Hazardous	-	-
RET/NTRK1 KIT RT-PCR Enzyme	Glycerol (Not Hazardous)	56-81-5	20-60%
RET/NTRK1 KIT Hy Buffer	Not Hazardous	-	-
SA-PE solution	Sodium azide	26628-22-8	0.05%
RET/NTRK1 KIT POS Control	Not Hazardous	-	-
NEG Control	Not Hazardous	-	-
	Sodium azide	26628-22-8	0.05%
RET/NTRK1 KIT πCode MicroDisc	Glycerol (Not Hazardous)	56-81-5	40-60%
10X Wash Buffer	Sodium azide	26628-22-8	0.09%
	Tween 20 (Not Hazardous)	9005-64-5	0.1%
ddH2O	Not Hazardous	-	-

Doc No.: B47-056-V03



Ver. 03 08/24/2020

IntelliPlex RET/NTRK1 Rearrangement Kit

The concentration or concentration ranges of all ingredients which are hazardous within the meaning of the GHS are present below their cut-off levels. No components need to be disclosed according to the applicable regulations.

Component	: Sodium azide	Classification
Formula	NaN ₃	Acute Tox. 2 (H300); Acute Tox. 1 (H310);
Cas-No	26628-22-8	STOT RE 2 (H373); Aquatic Acute 1 (H400);
EC-No	247-852-1	Aquatic Chronic 1 (H410)
Index-No.	011-004-00-7	(EUH032)
For the fu	Il text of the H-Statements me	entioned in this Section, see Section 16.

SECTION 4: FIRST AID MEASURES

4.1 Description of first aid measures		
General information:	Consult a physician. Show this safety data sheet to the doctor in attendance. Move out of dangerous area.	
Eye Contact:	Flush eyes with water as a precaution. Remove contact lenses. Consult a physician.	
Skin Contact:	Wash off with soap and plenty of water. Consult a physician.	
Inhalation:	If breathed in, move person into fresh air. If not breathing, give artificial respiration. Consult a physician	
Ingestion:	Never give anything by mouth to an unconscious person. Rinse mouth with water. Consult a physician.	
4.2 Most important syn	nptoms and effects, both acute and delayed	
The most important known symptoms and effects are described in section 2, section 3 and/or in section 11		
13 Indication of any im	mediate medical attention and special treatment needed	

4.3 Indication of any immediate medical attention and special treatment needed

No data available

SECTION 5: FIRE AND EXPLOSION DATA

5.1 Extinguishing media

Standard procedure for chemical fires. Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

- 5.2 Special hazards arising from the substance or mixture
- Not flammable or combustible.

5.3 Advice for firefighters

Wear self-contained breathing apparatus for firefighting if necessary.

5.4 Further information

Hazardous decomposition products formed under fire conditions: Sodium oxides.

SECTION 6: ACCIDENTAL RELEASE MEASURES

6.1 Personal precautions, protective equipment and emergency procedures

Wear respiratory protection. Avoid dust formation. Avoid breathing vapors, mist or gas. Ensure adequate ventilation. Evacuate personnel to safe areas. Avoid breathing dust. For handling and storage refer to section 7. For personal protection see section 8.

6.2 Environmental precautions

Prevent further leakage or spillage if safe to do so. Do not let product enter drains. Discharge into the environment must be avoided.

6.3 Methods and materials for containment and cleaning up

Take up mechanically and collect in suitable container for disposal.

6.4 Reference to other sections



Ver. 03 08/24/2020

For handling and storage refer to section 7. For personal protection see section 8. For disposal see section 13.

SECTION 7: HANDLING AND STORAGE

7.1 – Precautions for safe handling

Avoid contact with skin and eyes. Avoid formation of dust and aerosols. Provide appropriate exhaust ventilation at places where dust is formed.

7.2 – Conditions for safe storage, including anything that is incompatible

Keep container tightly closed in a dry and well-ventilated place. Refer to product label.

7.3 – Specific end use(s)

Refer to section 1.2

SECTION 8: EXPOSURE CONTROLS AND PERSONAL PROTECTION

8.1 Control parameters **United States of America**

Component	ACGIH TLV	OSHA PEL	NIOSH IDLH
Sodium azide	Ceiling: 0.29 mg/m ³ Ceiling: 0.11 ppm	Skin (Vacated) Ceiling: 0.1 ppm (Vacated) Ceiling: 0.3 mg/m³	Ceiling: 0.1 ppm Ceiling: 0.3 mg/m ³
Glycerol		(Vacated) TWA: 10 mg/m ³ (Vacated) TWA: 5 mg/m ³ TWA: 15 mg/m ³ TWA: 5 mg/m ³	
TS-769099 (RNAstable LD)	TWA: 10 ppm 25 mg/m³	TWA: 10 ppm 25 mg/m ³	

ACGIH - American Conference of Governmental Industrial Hygienists OSHA - Occupational Safety and Health Administration

NIOSH IDLH: NIOSH - National Institute for Occupational Safety and Health

European Union

Component	CAS-No.	Value type (Form	Control parameters	Basis
		of exposure)		
Sodium Azide	26628-22-8			
		TWA	0.1 mg/m3	2000/39/EC and 2006/15/EC
Further information	Identifies the	possibility of significan	t uptake through the sk	kin, Indicative
		STEL	0.3 mg/m3	2000/39/EC and 2006/15/EC
Further information	Identifies the possibility of significant uptake through the skin, Indicative			
		AGW	0.2 mg/m3	DE TRGS 900
Glycerol	56-81-5			
		AGW (Inhalable fraction)	200 mg/m3	DE TRGS 900
		TWA	10 mg/m3	TLV-ACGIH

8.2 Exposure controls

Appropriate engineering controls:



Ver. 03 08/24/2020

Use process enclosures, local exhaust ventilation, or other engineering controls to keep airborne levels below recommended exposure limits. If user operations generate dust, fume or mist, use ventilation to keep exposure to airborne contaminants below the exposure limit.

Personal protective equipment

Follow usual standard laboratory practices. Use appropriate chemical resistant gloves, appropriate safety glasses and wear protective work clothing.

Eye Protection

Wear appropriate protective eyeglasses or chemical safety goggles as described by OSHA's eye and face protection regulations in 29 CFR 1910.133 or European Standard EN166.

Hand Protection

Protective gloves (Natural rubber, Nitrile rubber, Neoprene, PVC; Breakthrough time see manufacturers information) **Skin and body protection**

Long sleeved clothing

Respiratory protection

When working with concentrations above the exposure limit use engineering controls and air-purifying respirators (full-face particle respirator type N100 (US) or type P3 (EN 143)). Use a full-face supplied air respirator if no engineered controls are utilized. Use respirators and components tested and approved under appropriate government standards such as NIOSH (US) or CEN (EU).

Environmental exposure controls

Prevent leakage or spillage. Discharge into the environment must be avoided.

SECTION 9: PHYSICAL AND CHEMICAL PROPERTIES

9.1 Information on basic physical and chemical properties

RET/NTRK1 KIT RT-PCR Reaction Mix; RET/NTRK1 KIT RT-PCR Primer; RET/NTRK1 KIT RT-PCR Hy Buffer; SA-PE solution; RET/NTRK1 KIT RT-PCR POS Control: NEG Control: RET/NTRK1 KIT RT-PCR πCode MicroDisc: 10X Wash Buffer

Physical State and Appearance:	Liquid.	Vapor Pressure:	Not available.
Odor:	Not available.	Vapor Density:	Not available.
Taste:	Not available.	Volatility:	Not available.
pH (1% soln/water):	Not available.	Odor Threshold:	Not available.
Boiling Point:	Not available.	Water/Oil Dist. Coeff.:	Not available.
Melting Point:	Not available.	Ionicity (in Water):	Not available.
Critical Temperature:	Not available.	Dispersion Properties:	Not available.
Specific Gravity:	Not available.	Solubility:	Not available.
Physical State and Appearance:	Solid.	Vapor Pressure:	Not available
	Solid. Not available.	Vapor Pressure: Vapor Density:	Not available Not available
Odor:		-	
Odor: Taste:	Not available.	Vapor Density:	Not available
Physical State and Appearance: Odor: Taste: pH (1% soln/water): Boiling Point:	Not available. Not available.	Vapor Density: Volatility:	Not available Not available
Odor: Taste: pH (1% soln/water): Boiling Point:	Not available. Not available. Not available.	Vapor Density: Volatility: Odor Threshold:	Not available Not available Not available
Odor: Taste: pH (1% soln/water):	Not available. Not available. Not available. Not available.	Vapor Density: Volatility: Odor Threshold: Water/Oil Dist. Coeff.:	Not available Not available Not available Not available

9.2 Other safety information

Not available



Ver. 03 08/24/2020

SECTION 10: STABILITY AND REACTIVITY DATA				
10.1 Reactivity:	The product is stable under recommended shipping and storage conditions.			
10.2 Chemical Stability:	The product is stable under recommended shipping and storage conditions.			
10.3 Possibility of hazardous reactions	Will not occur.			
10.4 Conditions to Avoid	Excess heat			
10.5 Incompatible materials 10.6 Hazardous decomposition products	Strong acids/alkalis, strong oxidizing/reducing agents. Not available.			

SECTION 11: TOXICOLOGICAL INFORMATION		
11.1 Information on toxicological effects		
Acute toxicity:	Not classified based on available information.	
Ingredients		
Sodium Azide	Acute toxicity LD50 Oral - Rat - 27 mg/kg	
Glycerol	Acute toxicity LD50 Oral - Rat - 12600 mg/kg	
	LD50 Dermal - Rabbit - > 10000 mg/kg	
Tween20	LD50 Oral - Rat - 36700 mg/kg	
Eye Contact:	No data available	
Skin Contact:	No data available	
Inhalation:	No data available	
Ingestion:	No data available	
Carcinogenic Effects:	IARC: No component of this product present at levels greater than or	
	equal to 0.1% is identified as probable, possible or confirmed human carcinogen by IARC.	
	ACGIH: No component of this product present at levels greater than or	
	equal to 0.1% is identified as a carcinogen or potential carcinogen by ACGIH.	
	NTP: No component of this product present at levels greater than or	
	equal to 0.1% is identified as a known or anticipated carcinogen by NTP.	
	OSHA: No component of this product present at levels greater than or	
	equal to 0.1% is on OSHA's list of regulated carcinogens.	
Mutagenic Effects:	No data available	
Reproduction Toxicity:	No data available	
Sensitization:	No data available	
STOT - Single exposure	No data available	
STOT - Repeated exposure	No data available	
Additional Information:	RTECS: VY8050000; MA8050000; TR7400000	

SECTION 12: ECOLOGICAL INFORMATION		
12.1 Toxicity:	No data available	
12.2 Persistence and degradability:	No data available	
12.3 Bioaccumulative potential:	No data available	
12.4 Mobility in soil:	No data available	
12.5 Results of PBT and vPvB assessment:	This substance/mixture contains no components considered to be either persistent, bioaccumulative and toxic (PBT), or very persistent	



SAFETY DATA SHEET

IntelliPlex RET/NTRK1 Rearrangement Kit

Ver. 03 08/24/2020

12.6 Other adverse effects:

and very bioaccumulative (vPvB) at levels of 0.1% or higher.

An environmental hazard cannot be excluded in the event of unprofessional handling or disposal. Potentially toxic to aquatic life.

SECTION 13: DISPOSAL CONSIDERATIONS

13.1 Waste treatment methods

Product: Offer surplus and non-recyclable solutions to a licensed disposal company. **Contaminated packaging:** Dispose of as unused product.

SECTION 14: TRANSPORT INFORMATION

General

The product is not covered by international regulations on the transport of dangerous goods (IMDG, IATA, DOT, ADR/RID).

14.1 – UN Number	Not regulated Not applicable
14.2 – Proper shipping name	Not regulated Not applicable
14.3 – Hazard class	Not regulated Not applicable
14.4 – Packing group	Not regulated Not applicable
14.5 – Environmental hazards	Not regulated Not applicable
14.6 – Special precautions for the user	Not regulated Not applicable
14.7 – Transport in bulk information	Not regulated Not applicable

SECTION 15: OTHER REGULATORY INFORMATION

European Union:

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture Contains no REACH substances with Annex XVII restrictions. Contains no substance on the REACH candidate list. Contains no REACH Annex XIV substances.

15.2 Chemical Safety Assessment

For this product a chemical safety assessment was not carried out

United States of America:

Regulatory Status Classified in accordance with Appendix A, Appendix B and Appendix F of the OSHA Hazard Communication Standard 29 CFR §1910.1200. Regulatory References OSHA Hazard Communication Standard 29 CFR §1910.1200

All components of this product are listed, or excluded from listing, on the United States Environmental Protection Agency Toxic Substances Control Act (TSCA) inventory

SARA 302 Components

No chemicals in this material are subject to the reporting requirements of SARA Title III, Section 302.

SARA 313 Components

This material does not contain any chemical components with known CAS numbers that exceed the threshold (De

Minimis) reporting levels established by SARA Title III, Section 313.

SARA 311/312 Hazards

See section 2 for more information

Clean Water Act/ Clean Air Act

This product does not contain any substances regulated as pollutants pursuant to the Clean Water Act or Clean Air Act.



Ver. 03 08/24/2020

US State Regulation

State Right to Know:

Glycerol (CAS-No. 56-81-5) can be found on the following state Right-to-Know lists: Pennsylvania, Minnesota, Massachusetts, New Jersey, Rhode Island

Sodium azide (CAS-No. 26628-22-8) can be found on the following state Right-to-Know lists: Pennsylvania, Minnesota, Massachusetts, New Jersey, Rhode Island

Tween 20 (Cas-No. 9005-64-5) can be found on the following state Right-to-Know lists: Pennsylvania, Minnesota, Massachusetts, New Jersey, Rhode Island

TS-769099 (RNAstable LB) can be found on the following state Right-to-Know lists: Pennsylvania, Minnesota, Massachusetts, New Jersey, Rhode Island

California Directors List of Hazardous Substances

The following ingredients are listed or exempt: TS-769099

California Prop. 65 Components

This product does not contain any chemicals known to State of California to cause cancer, birth defects, or any other reproductive harm.

SECTION 16: OTHER INFORMATION

Full text of H-Statements referred to under section 3.

H300 Fatal if swallowed.

H310 Fatal in contact with skin.

H373 May cause damage to organs through prolonged or repeated exposure if swallowed.

H400 Very toxic to aquatic life.

H410 Very toxic to aquatic life with long lasting effects.

EUH032 Contact with acids liberates very toxic gas.

To the best of our knowledge, the information contained herein is accurate and complete. However, we can neither guarantee nor assume any liability whatsoever for the accuracy or completeness of the information contained in this SDS. Final determination of suitability of any material is the sole responsibility of the user, as health and safety precautions in this data sheet may not be adequate for all individuals and/or situations. It is the user's obligation to evaluate and use this product safely and to comply with all applicable laws and regulations. All materials and mixtures may present unknown hazards and should be used with caution. No warranty is made, either express or implied.